

Medical Patient Weight Scale

Model : WM302

Medical patient weight scale with a maximum capacity of 300 kg. The WM302 has an extra large platform and handrails for unsteady patients. Ideal for medical practitioners, hospitals and retirement homes.

FEATURES

WM302

- professional patient scale ideal for hospitals, medical clinics and retirement homes
- sturdy, large diameter handrail design with rubber grip provides extra support and safety during patient weighing
- easy-to-read 25 mm LCD display with backlight
- auto weight-lock function on stable weight
- weight hold and last weight recall function
- body mass index (BMI) calculation upon height input
- wheels at rear of platform allow for easy relocation
- heavy duty low profile weigh platform with solid robust steel construction
- AC/DC adaptor and built in rechargeable battery operation
- battery life approximately 60 hours continuous use
- Australian Register of Therapeutic Goods (ARTG) # 174843

SPECIFICATIONS

Capacity:	300 kg x 100 g
Platform size (mm):	600 (W) x 600 (D) x 60 (H)
Dimensions (mm):	600 (W) x 650 (D) x 950 (H)
Net Weight:	23.4 kg

OPTIONS

- integrated height rod kit with measuring range of 65 – 211 cm x 1 mm increments
- WMTP printer



Specifications may change without notice