


**Wedderburn**<sup>®</sup>  
Since 1896.

## Patient Transfer Scale

Product Code: MDPTSM999

This durable patient transfer board with integrated weighing scale provides a revolutionary solution for weighing time-critical or immobile patients quickly, safely and accurately.



## Patient Transfer Scale

for weighing time-critical or immobile patients safely, quickly and accurately.

Weighing critical patients is a crucial process and also a challenging one.

Weighing patients correctly is essential to ensuring accurate drug dosage and administering the correct treatment. For patients requiring urgent immediate medical attention, accurate weighing is not always easily achievable. It's for this reason, and through the determination of an acute nurse working for NHS Lanarkshire to find a reliable solution, that the Patient Transfer Scale was developed.

Key benefits include:

- ▲ One device for lateral transfer and weighing of patients
- ▲ Weigh immobile patients in seconds with less stress and more dignity for the patient
- ▲ Accurate weight measurement when it's critical for safe medication management and patient monitoring



### Instant Accurate Weighing when you need it most!

When knowing the weight of a patient is absolutely crucial to delivering the right care and treatment, the MDPTSM999 embeds confidence and accuracy in your processes.

Use the PTS throughout your hospital pathway to enhance outcomes particularly for; trauma patients in emergency departments where time is crucial, stroke patients to ensure a quick 'door to needle' time, plus Cardiac, ICU, HTU, and X-ray patients.

### Time Saving Integrated Process

Our Patient Transfer Scale integrates a transfer board and patient weighing scale into one unit.

This means weighing patients becomes part of the already existing process of transferring a patient from a stretcher to bed, or bed to bed.

It's quick and easy, and saves you valuable time.



# Patient Transfer Scale

MDPTSM999



## Easy to Read Display

The bright LCD is easy to read and can be programmed to automatically power-off after 30 or 60 minutes.

The backlight can easily be set to ON or OFF making the display even easier to read whatever lighting conditions you are in.

## Lightweight and Portable

The MDPTSM999 has been made as lightweight as possible so it's easy to take it to exactly where you need it and get an instant, accurate weight.


Grab handles are conveniently positioned down each side making it easier to move.



## Features Overview

- lightweight, durable construction
  - low profile (30 mm)
- easy to clean wipe down surface
- 10 conveniently placed grab handles
- powered by rechargeable battery
  - weight hold function
  - bright easy to read LCD display
- display options include programmable auto-off and backlight
- easy to use so requires minimal training

## Features Overview



### SPECIFICATIONS

ARTG Number	315361
Maximum Capacity	250 kg
Graduation	0.5 kg
Weighing Units	kg
Type of Display	27.7 mm x 75 mm LCD
Environment Suitability	Dry
Function Keys	ON/OFF/ZERO, HOLD
Operating Temperature	5 to 35 °C
Power Source	Rechargeable battery pack (up to 50 hours continuous use )
Net Weight	11 kg
Weighing Area (mm)	1800 (W) x 510 (D)
Overall Dimensions (mm)	1800 (W) x 700 (D) x 30 (H)

*Specifications may change without notice*