

Medical Patient Weight Scale

Model : WM202

Medical patient weight scale with a maximum weighing capacity of 250 kg.
Compact and light weight non-slip platform with solid steel structure.
Ideal for medical practitioners throughout Australia.

FEATURES

WM202

- professional patient scale ideal for hospitals, medical clinics and health clubs
- easy-to-read 25 mm LCD display with backlight
- auto weight-lock function on stable weight
- weight hold and last weight recall function
- Body Mass Index (BMI) calculation upon height input
- compact light weight ABS non slip platform, with solid steel structure
- AC/DC adaptor with built in rechargeable battery operation
- battery life approximately 60 hours continuous use
- net weight 10.5 kg
- Australian Register of Therapeutic Goods (ARTG) # 174843

OPTION

- height rod kit with measuring range of 65 – 211cm x 1mm increment



Shown with optional height rod

SPECIFICATIONS

Model Number	Max Capacity	Graduation	Platform Size (mm)	Overall Dimensions (mm)
WM202	250 kg	100 g	340 (W) x 440 (D) x 120 (H)	340 (W) x 440 (D) x 1030 (H)

Specifications may change without notice

▲ **Weighing Scales.**

*Labels, Labelling Equipment, Packaging Equipment, Weighing Scales,
POS Systems, Food Equipment, Maintenance & Repairs, Hire.*