

Medical Patient Chair Scale

Model : WM403

Medical patient chair scale with 300 kg capacity, retractable armrests and foldable footrest with an easy to read LCD display with backlight. Ideal for hospitals, retirement homes and medical clinics.

FEATURES

WM403

- high capacity professional patient chair scale ideal for hospitals, medical clinics and retirement homes
- sturdy steel tubular frame; powder coated for ease of cleaning
- retractable armrests and foldable footrest allowing unimpeded chair access
- moulded ergonomically designed ABS plus size seat providing comfort while easily accommodating larger patients
- effortless manoeuvrability with front and rear swivel wheels, rear wheels lockable
- easy to read 25 mm LCD display with backlight
- auto weight-lock function on stable weight
- weight hold and last weight recall function
- Body Mass Index (BMI) calculation upon height input
- AC/DC adaptor and battery operation (6 x AA batteries, batteries not included)
- battery life approximately 75 hours continuous use (alkaline) without backlight
- net weight 23.5 kg
- Australian Register of Therapeutic Goods (ARTG) # 174844



SPECIFICATIONS

| Model Number | Max Capacity | Graduation | Seat Size (mm) | Overall Dimensions (mm) |
|--------------|--------------|------------|-------------------|-----------------------------|
| WM403 | 300 kg | 100 g | 500 (W) x 480 (D) | 565 (W) x 840 (D) x 940 (H) |

Specifications may change without notice

▲ **Weighing Scales.**

Labels, Labelling Equipment, Packaging Equipment, Weighing Scales, POS Systems, Food Equipment, Maintenance & Repairs, Hire.