

# Medical Patient Wheelchair Scale

Model : WM501

Medical patient wheelchair scale. The WM501 has a low profile platform design which makes it easy to roll on and off when using a wheelchair. The platform is large enough to handle most types of wheelchairs. This multi-purpose scale is ideal for hospitals, medical clinics and retirement homes.

## FEATURES

WM501

- professional patient, multi purpose scale ideal for hospitals, medical clinics and retirement homes
- heavy duty low profile steel construction with integrated ramps affording easy on-off access for patients
- sturdy rear mounted hand rail and integrated ramps ensures patient safety when weighing either standing, with a walker, in wheel chairs or on a stationary chair
- collapsible hand rail, two lift handles and wheels at the rear of weigh platform allow for easy relocation and / or space saving storage
- easy to read 25 mm LCD display with back light
- auto weight-lock function on stable weight
- weight hold and last weight recall function
- Body Mass Index (BMI) calculation upon height input
- AC/DC adaptor and built in rechargeable battery operation
- battery life approximately 60 hours continuous use
- net weight 37 kg
- Australian Register of Therapeutic Goods (ARTG) # 174863



optional expanding ramps



## OPTIONS

- auxilliary expanding ramps allowing smooth ease of access to the weigh platform
- optional WMTTP printer

## SPECIFICATIONS

Model Number	Max Capacity	Graduation	Platform Size (mm)	Overall Dimensions (mm)
WM501	350 kg	100 g	1080 (W) x 990 (D) x 60 (H) Flat Area: 760 (W) x 800 (D)	1080 (W) x 990 (D) x 1000 (H)

Specifications may change without notice