

Wheelchair Scale - Portable

Model : WM502

Portable wheelchair scale with 300 kg capacity and BMI function. Integrated ramps on both ends allows for easy on and off access. This scale is ideal for hospitals, retirement homes and medical clinics where portability is required or where there is limited space.

FEATURES

WM502

- compatible with most wheelchairs, easily transported and ideal for home visits
- low angle ramps provide easy on-off access
- indicator can be placed on a desk or wall mounted
- easy to read 25 mm high LCD digital display
- improves weighing accuracy for unsteady or shaking patients
- tare function to subtract the weight of wheelchairs and patient's belongings
- improves record keeping by displaying the patient weight after they have left the scale
- Body Mass Index (BMI) function
- AC power adaptor or alkaline battery (not included) powered
- net weight 19 kg
- Australian Register of Therapeutic Goods (ARTG) #174863



Low angle ramps provide easy on-off access



Image shows optional carry case

OPTION

- carry case



SPECIFICATIONS

Model Number	Max Capacity	Graduation	Live Weighing Area (mm)	Overall Beam Dimensions (mm)*
WM502	300 kg	100 g	920 (L) x 130 (W)	1180 (L) x 170 (W) x 75 (H)

* Note, beam dimensions are for a single beam and include the ramp and handle

Specifications may change without notice