

Professional Personal Scale

Model : TIWB380H

The WB380H is a versatile, low maintenance precision personal scale with a high 300kg capacity and built in height rod. Its large low profile platform makes it easy to step on the scale, along with the built-in BMI function and rotating display making the WB380H ideal for use in a variety of applications.

FEATURES

TIWB380H

- high 300kg capacity
- built in height rod 61cm-250cm with 1cm increments
- large easy to read LCD display
- battery or AC/DC adaptor power
- large and ultra-low profile platform
- BMI function
- step-on activation
- auto power off
- up to 100hrs of continuous use battery life

SPECIFICATIONS

Dimensions (mm):	395(W) x 552(D) x 1425(H)
Platform dimensions (mm):	395(W) x 390(D) x 59(H)
Capacity:	300kg x 100g
Net weight (kg):	11.2

OPTIONS

- GMON Software



Specifications may change without notice