

# Professional Personal Scale

Model : TIWB380NP

The TIWB380NP is a low maintenance personal scale that delivers precision accuracy whilst weighing up to a significant 300kg capacity. It is low maintenance and very easy to use and features a large, stable yet very low-profile platform which makes it easy and safe for your clients and patients to step on. This highly accurate personal scale has a remote display, built-in BMI functionality and is very portable which makes it incredibly versatile and ideal for use in a variety of applications.

## FEATURES

TIWB380NP

- high 300kg capacity
- large easy to read LCD display
- battery or AC/DC adaptor power
- large and ultra-low profile platform
- BMI function
- step-on activation
- auto power off
- up to 100hrs of continuous use battery life
- display cable length 1 meter
- USB and RS232 interface

## SPECIFICATIONS

Platform size (mm):	395 (W) x 390 (D) x 59 (H)
Indicator size (mm):	240 (W) x 123 (D) x 139 (H)
Capacity:	300kg x 0.1kg
Net Weight:	7.1kg

## OPTIONS

- GMON software
- available pole mounted display model (TIWB380P)



*Specifications may change without notice*