

Professional Personal Scale

Model : TIWB380P

The WB380P is a versatile, low maintenance precision personal scale with a high 300kg capacity. Its large low profile platform makes it easy to step on the scale, along with the built-in BMI function and rotating display making the WB380P ideal for use in a variety of applications.

FEATURES

TIWB380P

- high 300kg capacity
- large easy to read LCD display
- battery or AC/DC adaptor power
- large and ultra-low profile platform
- BMI function
- step-on activation
- auto power off
- up to 100hrs of continuous use battery life
- USB and RS232 interface

SPECIFICATIONS

Platform size (mm):	395 (W) x 390 (D) x 59 (H)
Overall size (mm):	395 (W) x 595 (D) x 1196 (H)
Capacity:	300kg x 0.1kg
Net Weight:	10 kg

OPTIONS

- GMON software
- available model without pole (TIWB380NP)



Specifications may change without notice