

Height Rod with Stand

Model : WS220S

The Wedderburn WS-220S Height Rod with Stand is a free standing, stadiometer allowing a standing height measurement to be recorded with an accuracy of 1 mm.

From its high quality stainless steel base plate and stand, with aluminium telescopic measuring upright and durable ABS measuring head rest it has been designed for many years service.

With a maximum measuring range of 200cm, and its simplicity of use makes it ideal for the height measurement of children through to adults.

Wedderburn has long been the professional's choice for height measurement, due to our long trusted reliability and accuracy.



SPECIFICATIONS

WS220S

Measuring Range:	80 cm - 200 cm
Graduation:	1 mm
Net Weight:	4 kg

Specifications may change without notice