

Baby Scale

Model : WMMS3500

The WMMS3500 infant/baby scale features a curved baby tray to safely and comfortably hold infants during weighing. With 5 gram graduations below 10 kg the scale is able to track even slight changes in infants' weight. With a removable baby tray and 20 kg capacity it can also be used to weigh toddlers.

FEATURES

WMMS3500

- large easy to read LCD display
- simple three button operation
- large ergonomically contoured baby tray for added safety
- removable baby tray for weighing toddlers
- AC/DC adaptor or 9 V battery operation
- net weight 2.6 kg
- Australian Register of Therapeutic Goods (ARTG) # 329077

OPTIONS

- baby measuring rod, model # WMHM80M3500 (ARTG # 329078)
- scale carry bag as pictured, model # WSBACKPACK (not TGA registered)



Optional carry bag



VARIATIONS

Model Number	Max Capacity	Graduation	Baby Tray (mm)	Overall Dimensions (mm)
WMMS3500	20 kg	0 - 10 kg x 5 g 10 kg - 20 kg x 10 g	560 (W) x 290 (D) x 65 (H)	560 (W) x 350 (D) x 120 (H)

Specifications may change without notice