

Patient Scale

Model : WMMS6111

The WMMS6111 adult and child scale delivers outstanding performance and precision both for adult and child patients.

The low profile platform makes it easy to step on and off the scale, and a dual-sided LCD display allows the weight to be read by the medical professionals or the patient.

FEATURES

WMMS6111

- Body Mass Index (BMI) function
- hold function
- tare function
- reversible display
- combined infant weighing mode (10 gram graduation under 20 kg) and adult weighing mode (100 gram graduation between 20 - 250 kg)
- AC/DC adaptor and battery operation (4 x AA size batteries, batteries not included)
- suitable for use on flat hard surfaces
- net weight 2.5 kg
- Australian Register of Therapeutic Goods (ARTG) # 329076



SPECIFICATIONS

Model Number	Max Capacity	Graduation	Platform Size (mm)	Overall Dimensions (mm)
WMMS6111	250 kg	0 - 20 kg x 10 g 20 kg - 250 kg x 100 g	327 (W) x 260 (D)	345 (W) x 350 (D) x 45 (H)

Specifications may change without notice

▲ Weighing Scales.

Labels, Labelling Equipment, Packaging Equipment, Weighing Scales,
POS Systems, Food Equipment, Maintenance & Repairs, Hire.