

BC-418



NULL MODEM



graues Kabel



grey cable

PC-COM (RS 232)



MC-180, WB-100/110 BC-420, SC-330



ViScan



graues Kabel



grey cable



BF-350




TBF-300 ... 410



rotes Kabel



red cable




MC-180 MC-980 MC-780 SC-240



USB - Kabel
USB - cable



USB / RS232
Adapter



USB-Gerätetreiber auf
GMON-CD im
Verzeichnis "driver"



USB-Driver on GMON-CD
at directory "device driver"