

Baby Measuring Mat

Model : WMHM110M

Measuring baby length is as easy as laying them down on top of the HM110M Infant Measuring Mat.

The soft, flexible, vinyl mat is printed on both edges with measurements available in centimetres and inches. The rigid footboard folds to the bottom of the infant's feet and works with the headboard to ensure accurate measurements, with instructions printed directly on the mat.

FEATURES

WMHM110M

- graduation in centimetres and inches along both edges for easy viewing
- rounded corners for baby's safety
- non-toxic and safe for baby's skin
- can be cleaned with medical wipes
- high-quality screen print
- instructional diagram printed directly on mat
- Australian Register of Therapeutic Goods (ARTG) # 329078

SPECIFICATIONS

Measurement range:	10 cm - 110 cm x 5 mm
Overall dimensions (mm):	1400 (L) x 300 (W)
Weight:	0.4 kg



Specifications may change without notice