

Neonatal Baby Scale

Model : WM101

Ergonomically designed and fully featured Neonatal baby scale. Features include no tilt design, high side walls for added safety, auto hold weighing function, last weight recall, memory for up to nine individual weights, weight gain feature and hygienic ABS plastic construction. Ideal for use in Hospitals and Neonatal Care facilities.

FEATURES

WM101

- stable no tilt design with high tray walls for safe weighing
- large hygienic and easy to clean platter
- easy to read 23 mm high LCD display
- auto weight hold function provides accuracy even with moving babies
- weight gain function
- last weight recall function
- memory for up to 9 individual weights
- tare function to subtract weight of blanket or nappies
- low battery indicator
- AC/DC adaptor or 6 x AA battery operation
- net weight 4.6 kg
- Australian Register of Therapeutic Goods (ARTG) # 190276

OPTION

- scale carry bag



SPECIFICATIONS

Model Number	Max Capacity	Graduation	Baby Tray Dimensions (mm)	Overall Dimensions (mm)
WM101	15 kg	0 - 6 kg x 2 g 6 kg - 16 kg x 5 g	600 (W) x 270 (D) x 130 (H)	600 (W) x 300 (D) x 160 (H)

Specifications may change without notice