

Tanita Health Monitor Software

Model : TIGMONSW

Easily connected via RS232C or USB, the Tanita Health Monitoring software is a database capable of storing client information together with their measurements. It will also provide professional reports and printouts, trend analysis of composition measurements against time and much more. The software is an essential tool for use in health care, weight management, and sports and leisure facilities.



FEATURES

TIGMONSW

- designed for use with Tanita's professional range of body compositional equipment
- client data is automatically transferred from PC to body composition monitor
- generate a variety of reports and printouts
- chart client progress and trends
- ideal for use in health screening, weight management and fitness assessment
- connected via RS232C or USB
- compatible with Windows operating systems

SUITABLE FOR USE WITH TANITA MONITORS

- | | | |
|--------|----------|-----------|
| SC-240 | BC-418 | WB-110AZS |
| SC-330 | MC-780 | DC-360 |
| BF-410 | MC-980 | DC-430 |
| BC-420 | WB-3000P | WB-380 |



*Computer not included with software

Specifications may change without notice

▲ Weighing Scales.

Labels, Labelling Equipment, Packaging Equipment, Weighing Scales, POS Systems, Food Equipment, Maintenance & Repairs, Hire.



Provide professional reports and printouts, trend analysis of composition measurements against time and much more.



Specifications may change without notice

▲ Weighing Scales.

Labels, Labelling Equipment, Packaging Equipment, Weighing Scales, POS Systems, Food Equipment, Maintenance & Repairs, Hire.